

Enjoy deep restorative sleep with a soft cloth mask.

“Another hit! It is like you are wearing a pillow of air. A must try.”

Brian S.
CPAP User
Lawrence, KS

“Better than anything I have used!”

Larry H.
CPAP User
Gravette, AR

“I really love this mask! So much more comfortable than my old plastic mask.”

Theresa I.
CPAP User
Phoenix, AZ

“Best mask and headgear on the market.”

Barbara H.
CPAP User
Monroe Township, NJ

“Before SleepWeaver, I had given up on the CPAP machine.”

Susan L.
CPAP User
Astor, FL



Open Face Head-Strap Design:

Allows virtually unrestricted vision; especially good for people who like to fall asleep while reading or watching TV, accommodates eyeglasses.

90° Swivel Connector with 360° Rotation:

Allows users to sleep in more natural positions; easy to remove and clean.

Bottom Exhalation:

Exhalation holes are underneath the mask, approximating natural nasal exhalation. Angled so it doesn't bother users or bed partners.

Feather Weight Tube:

Lightweight 18" extension tube adds convenient quick release from standard hose; reduces drag on the soft cloth mask.

Super Lightweight:

Users barely notice the weight of this mask.

Sizes & Styles:

A number of different size options allows users to find the perfect fit for their face as well as their personality; expands the number of users the mask will fit.

Moisture Vapor Breathable:

High-tech material allows the skin to breathe naturally, wicks away moisture and leaves no pressure marks that are common with rubber or silicone seals.

Loose Fit:


As with all SleepWeaver masks, the unique “loose fit” provides the ultimate comfort and eliminates strap marks.

Tether Strap:

Allows users to position the tube over their head with the option to route the tube over the bed's headboard.

Easy Mask / Cushion Replacement:

Quick, simple and cost-effective replacement of material that routinely comes in contact with the skin surface.


















Available Colors




Starter Kit Colors

SleepWeaver Élan Ordering Information

Products	Photos	Sizes	HCPCS Codes	Part Numbers	Color Options			
SleepWeaver Élan Starter Kit		Large Regular Small	A7034 - Mask A7035 - Headgear A7037 - Tube	Starter Kit	100563	Beige		
				Starter Kit	100782	Blue		
SleepWeaver Élan Mask and Headgear		Large Regular Small	A7034 - Mask A7035 - Headgear	Large	Regular	Small		
				100481	100518	100727	Beige	
				100482	100519	100728	Blue	
SleepWeaver Élan Mask/Cushion		Large Regular Small	A7034 - Mask	Large	Regular	Small		
				100663	100653	100742	Beige	
				100664	100654	100743	Blue	
				100665	100655	100744	Pink	
				100666	100656	100745	Camouflage	
100667	100657	100746	Leopard					
SleepWeaver Élan Headgear Only		Universal	A7035 - Headgear	Headgear	100630	Black		
SleepWeaver Feather Weight Tube		18 inches	A7037 - Tube	Tube	100477	White		

SleepWeaver Élan Accessories

			
Tether Strap - 100638	Hose Clip - 100775	Connector - 100777	Washer - 100778
Tether Strap attaches to Tube & Headgear	Hose Clip prevents Tube movement	SleepWeaver Élan Threaded Connector Replacement Part	SleepWeaver Élan Washer Replacement Part

Contact us for more information about the SleepWeaver line of Soft Cloth PAP Masks:
 1-888-825-9640 * 1-724-858-2837
 orders@circadiance.com
 www.Circadiance.com



SleepWeaver
Soft Cloth PAP Masks

